





FLOOR

mm
MANACA
MEDICAL



LAUTUS
Antibacterial Mat

REF. 4571143MAT

LAUTUS ANTIBACTERIAL MAT

The multi-layer mat with an adhesive surface is effective for the control of contamination, bacteria and pathogens from the soles of footwear and wheels of transport carts. Helps reduce contamination before it enters clean rooms or critical environments.

FEATURES

- ✓ PE LD (low density polyethylene) multisheets of 0.05mm descending numbered with an 8µm thick acrylic adhesive with iodinated antibacterial/antifungal additive.
- ✓ Totally smooth numbered mat to prevent tripping hazard available in 30 layers with 10 mats per box.
- ✓ It is done by a consistent adhesive coverage ensuring maximum performance for each layer.
- ✓ Antimicrobial, resistant to mold and fungus.
- ✓ Automatically cleans 82% to 99.8% of bacteria from shoe soles, bed wheels, trolleys and stretchers that access clean rooms.





FLOOR

mm
MANACA
MEDICAL

REF. 4571143MAT

LAUTUS ANTIBACTERIAL MAT

SPECS

- ✓ Adhesiveness between sheets: 250 g / 25 mm \pm 40
- ✓ Heat resistance: 70° for 48 hours
- ✓ Adhesiveness to the ground: 500 g / 25 mm \pm 40
- ✓ Tensile strength 1.2 kg / cm

PRESENTATION

Measurements: 457 x 1143 mm

Block / Sheets: 30 layers with 10 mats / box

Transport box: 50 units. Each transport box includes 5 inner boxes of 10 units.



SINGLE USE



NON STERILE