



NON-WOVEN

mm
MANACA
MEDICAL



REF. KIT 2

BEDSHEETS SET

Set of single-use bed sheets that includes a 110x230 adjustable fitted sheet with elastic, a 260 x 160 top sheet and a 40 x 110 pillow cover.

FEATURES

- ✔ Disposable bed sheets, single use only.
- ✔ 25 grams polypropylene.
- ✔ Latex free.
- ✔ Hypoallergenic, non-irritating to skin and eyes.
- ✔ Air permeable.
- ✔ Hydrophobic (also blood repellent).
- ✔ Non-toxic and non-polluting when recycled.
- ✔ Dust-free.
- ✔ Low flammability (flame retardant).
- ✔ Heat resistant.
- ✔ Suitable for food contact.






NON-WOVEN

mm
MANACA
MEDICAL

REF. KIT 2

BEDSHEETS SET

Fitted sheet with elastic



Top sheet



PRESENTATION

- Fitted sheet with elastic 110 x 230 cm
- Top sheet 260 x 160 cm
- Pillow cover 40 x 110 cm



SINGLE USE



LATEX FREE



LOW FLAMMABILITY