



## STERILIZATION



REF. CTRV 100

### STERILIZATION COVERAGE

The Manaca Medical sterilization cover provides a wet and dry microbiological barrier according to the European directive EN868 part 1. And part 2 for non-woven fabric. It is not alcohol repellent so it does not meet the requirements to be used in the sterile field.

#### FEATURES

- ✓ External identification by adhesive labels on opposite corners.
- ✓ Polypropylene (SMS) formed by a tri-layer combination of Spunbond (S) long micro-metric fibers and Meltblown (M) short nanometric fibers.
- ✓ The Spunbond layers guarantee great flexibility and ease of bending, while the Meltblown layers provide an effective barrier against germs.
- ✓ Manufactured in 53 g/m<sup>2</sup> in purple color.

#### PRESENTATION

Box of 200 units, contained in a LD PET bag.  
Measurement of 100 x 100 cm.

#### COMPLIANCE

European directive EN868 part 1.



LATEX  
FREE



HEAT  
RESISTANT

