



SURGICAL

mm
MANACA
MEDICAL



REF. C8400(XS), C8401(S), C8302 (M), C8303 (L), C8304 (XL), C8305 (XXL)

WOUND RETRACTORS

Device that allows the surgeon to access the abdominal cavity through an atraumatically retracted wound, providing maximum exposure with minimum incision size. The Manaca Medical Wound Retractor is indicated for endoscopic and open surgery, especially for laparoscopic surgery, during removal of specimens in the absence of pneumoperitoneum. It provides clear exposure, protects wounds from contamination and reduces complications.

FEATURES

- ✓ The Manaca Medical disposable wound retractor is composed of an outer ring, an inner ring and a channel.
- ✓ Each part can be deformed longitudinally and radially by applying light force, upon cessation of force elastic recovery occurs.
- ✓ It provides maximum exposure for a clean operative field, the channel isolates the incision site from the abdominal viscera to reduce complications.
- ✓ Surgeons should use different wound retractors, depending on the site and size of the incision.





SURGICAL

REF. C8400(XS), C8401(S), C8302 (M), C8303 (L), C8304 (XL), C8305 (XXL)

WOUND RETRACTORS

REF.	SYMBOL	DESCRIPTION	INCISION SIZE	WOUND RETRACTION	THORACIC CAVITY RETRACTION	PNEUMOPERITONEUM MAINTENANCE	TROCAR ENTRY PORT
C8400	XS	Extra Small	2 to 4 cm	✓	✓	✓	✓
C8401	S	Small	2.5 to 6 cm	✓	✓	✓	✓
C8302	M	Medium	5 to 9 cm	✓	✓	✓	✓
C8303	L	Large	9 to 14 cm	✓			
C8304	XL	Extra Large	11 to 17 cm	✓			
C8305	XXL	Double Extra Large	17 to 25 cm	✓			

PRESENTATION

Individual box and individual sterile blister.

COMPLIANCE

Product in accordance with European standards of Regulation (EU) 2017/745, CE marked, ISO 13845,. Class II A sterile. Latex free and phthalate free.



SINGLE USE



LATEX FREE



PHthalate FREE



CE MARKING



STERILE